

HENLEY SELFRETAINING RETRACTOR ONLY

sharp; 3X4 teeth

16.0cm (6.4")




SKU: 070-300-16

Categories: [Open Surgery](#), [Retractors](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Selfretaining Retractors |
| AUTHORS | HENLEY |
| SPECIFICATIONS | sharp , toothed |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Selfretaining Retractors |
| AUTHORS | HENLEY |
| SPECIFICATIONS | sharp , toothed |