

HENLEY SELFRETAINING RETRACTOR ONLY

sharp; 3X4 teeth

16.0cm (6.4")




SKU: 070-300-16

Categories: [Open Surgery](#), [Retractors](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
AUTHORS	HENLEY
SPECIFICATIONS	sharp , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
AUTHORS	HENLEY
SPECIFICATIONS	sharp , toothed