

## HENLEY SELFRETAINING RETRACTOR

sharp; 3X4 teeth; complete with blades

16.5cm (6.6")




**SKU:** 070-300-00

**Categories:** [Open Surgery](#), [Retractors](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Selfretaining Retractors</a>
<b>KLS Martin</b>	15-769-00
<b>AUTHORS</b>	<a href="#">HENLEY</a>
<b>SPECIFICATIONS</b>	<a href="#">sharp</a> , <a href="#">toothed</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Selfretaining Retractors</a>
<b>KLS Martin</b>	15-769-00
<b>AUTHORS</b>	<a href="#">HENLEY</a>
<b>SPECIFICATIONS</b>	<a href="#">sharp</a> , <a href="#">toothed</a>