

HF BALL ELECTRODE; BIPOLAR; STERILE

Ø 2mm, single use

for 18.5 Fr. system; green; Box=5 pieces



SKU: 582-180-43B-S

Categories: [Endoscopy](#), [GUBBINI/Gyn](#), [Gynecology](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	Box
RANGE	Mini Resectoscopy (Midi)
TYPE	Electrodes

ZUSÄTZLICHE INFORMATIONEN

UNIT	Box
RANGE	Mini Resectoscopy (Midi)
TYPE	Electrodes