

## HF BALL ELECTRODE; BIPOLAR; STERILE

Ø 2mm, single use

for 18.5 Fr. system; green; Box=5 pieces



**SKU:** 582-180-43B-S

**Categories:** [Endoscopy](#), [GUBBINI/Gyn](#), [Gynecology](#), [Urology](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">Box</a>
<b>RANGE</b>	<a href="#">Mini Resectoscopy (Midi)</a>
<b>TYPE</b>	<a href="#">Electrodes</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">Box</a>
<b>RANGE</b>	<a href="#">Mini Resectoscopy (Midi)</a>
<b>TYPE</b>	<a href="#">Electrodes</a>