

HF BALL ELECTRODE; BIPOLAR; STERILE

Ø 3mm; single use; yellow/ blue

for 24/26 Fr. System; Box = 5 pieces



SKU: 630-240-15B-S

Categories: [Endoscopy](#), [Gynecology](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	Box
RANGE	Resectoscopy 4.0
TYPE	Electrodes

ZUSÄTZLICHE INFORMATIONEN

UNIT	Box
RANGE	Resectoscopy 4.0
TYPE	Electrodes