

## HF BALL ELECTRODE; BIPOLAR; STERILE

Ø 3mm; single use; yellow/ blue

for 24/26 Fr. System; Box = 5 pieces



**SKU:** 630-240-15B-S

**Categories:** [Endoscopy](#), [Gynecology](#), [Urology](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">Box</a>
<b>RANGE</b>	<a href="#">Resectoscopy 4.0</a>
<b>TYPE</b>	<a href="#">Electrodes</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">Box</a>
<b>RANGE</b>	<a href="#">Resectoscopy 4.0</a>
<b>TYPE</b>	<a href="#">Electrodes</a>