

HF-GENERATOR FOR ELECTRICAL SURGERY

for monopolar and bipolar frequency

PWT-400; Touch Screen

for bipolar resection



SKU: 540-280-50

Categories: [Endoscopy](#), [Equipment](#), [GUBBINI/Gyn](#), [Gynecology](#), [Laparoscopy](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [HF-Generator](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [HF-Generator](#)