

## HF KNIFE ELECTRODE; BIPOLAR; STERILE

angled; single use

for 16 Fr. System; black; Box = 5 pieces



**SKU:** 582-150-43GUB-S

**Categories:** [Endoscopy](#), [GUBBINI/Gyn](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">Box</a>
<b>RANGE</b>	<a href="#">Ellipse</a> , <a href="#">Original</a>
<b>TYPE</b>	<a href="#">Electrodes</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">Box</a>
<b>RANGE</b>	<a href="#">Ellipse</a> , <a href="#">Original</a>
<b>TYPE</b>	<a href="#">Electrodes</a>