

HF LOOP ELECTRODE; BIPOLAR; STERILE

90° angled; teardrop; single use

for 16 Fr. System; green; Box = 5 pieces



SKU: 582-150-40GUBE-S

Categories: [Endoscopy](#), [GUBBINI/Gyn](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	Box
RANGE	Ellipse
TYPE	Electrodes

ZUSÄTZLICHE INFORMATIONEN

UNIT	Box
RANGE	Ellipse
TYPE	Electrodes