

HIGH PERFORMANCE LIGHT GUIDE CABLE

for LED light source; w/o connectors

Ø 3.5 x 1800 mm Length

tube colour lilac

HP
HIGH PERFORMANCE LINE



SKU: 540-823-18LL

Categories: [Arthroscopy](#), [Endoscopy](#), [Equipment](#), [GUBBINI/Gyn](#), [Gynecology](#), [Laparoscopy](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Light Guide Cables
TYPE	HP (High Performance)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Light Guide Cables
TYPE	HP (High Performance)