

## HIGH PERFORMANCE LIGHT GUIDE CABLE

for LED light source; w/o connectors

Ø 3.5 x 1800 mm Length

tube colour yellow

**HP**  
HIGH PERFORMANCE LINE



**SKU:** 540-823-18YW

**Categories:** [Arthroscopy](#), [Endoscopy](#), [Equipment](#), [GUBBINI/Gyn](#), [Gynecology](#), [Laparoscopy](#), [Urology](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Light Guide Cables](#)

**TYPE** [HP \(High Performance\)](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Light Guide Cables](#)

**TYPE** [HP \(High Performance\)](#)