


HIGH PERFORMANCE LIGHT GUIDE CABLE

for LED light source; w/o connectors

Ø 3.5 x 2300 mm Length

tube colour lilac

HP
HIGH PERFORMANCE LINE



SKU: 540-823-23LL

Categories: [Arthroscopy](#), [Endoscopy](#), [Equipment](#), [GUBBINI/Gyn](#), [Gynecology](#), [Laparoscopy](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|--------------|---------------------------------------|
| UNIT | pcs. |
| RANGE | Light Guide Cables |
| TYPE | HP (High Performance) |

ZUSÄTZLICHE INFORMATIONEN

| | |
|--------------|---------------------------------------|
| UNIT | pcs. |
| RANGE | Light Guide Cables |
| TYPE | HP (High Performance) |