

HIGH PERFORMANCE LIGHT GUIDE CABLE

for LED light source; w/o connectors

Ø 3.5 x 2300 mm Length

tube colour yellow

HP
HIGH PERFORMANCE LINE



SKU: 540-823-23YW

Categories: [Arthroscopy](#), [Endoscopy](#), [Equipment](#), [GUBBINI/Gyn](#), [Gynecology](#), [Laparoscopy](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Light Guide Cables](#)

TYPE [HP \(High Performance\)](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Light Guide Cables](#)

TYPE [HP \(High Performance\)](#)