

## HIGH PERFORMANCE LIGHT GUIDE CABLE

for LED light source; w/o connectors

Ø 2.5 x 1800 mm Length

tube colour yellow

**HP**  
HIGH PERFORMANCE LINE



**SKU:** 540-822-18YW

**Categories:** [Arthroscopy](#), [Endoscopy](#), [Equipment](#), [GUBBINI/Gyn](#), [Gynecology](#), [Laparoscopy](#), [Urology](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Light Guide Cables</a>
<b>TYPE</b>	<a href="#">HP (High Performance)</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Light Guide Cables</a>
<b>TYPE</b>	<a href="#">HP (High Performance)</a>