

## HOFFMANN TUBING CLAMP

with fixation screw

25mm



**SKU:** 250-404-25

**Categories:** [Open Surgery](#), [Various Items](#)

### ZUSÄTZLICHE INFORMATIONEN

|                   |                               |
|-------------------|-------------------------------|
| <b>UNIT</b>       | <a href="#">pcs.</a>          |
| <b>TYPE</b>       | <a href="#">Tubing Clamps</a> |
| <b>KLS Martin</b> | 63-424-25                     |
| <b>AUTHORS</b>    | <a href="#">HOFFMANN</a>      |

### ZUSÄTZLICHE INFORMATIONEN

|                   |                               |
|-------------------|-------------------------------|
| <b>UNIT</b>       | <a href="#">pcs.</a>          |
| <b>TYPE</b>       | <a href="#">Tubing Clamps</a> |
| <b>KLS Martin</b> | 63-424-25                     |
| <b>AUTHORS</b>    | <a href="#">HOFFMANN</a>      |