

HOHMANN BONE LEVER

8mm; narrow tip, 1 whole handle

22.0cm (8.8")


SKU: 190-572-08

Categories: [Open Surgery](#), [Orthopaedics](#)



ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|-----------------------------|
| UNIT | pcs. |
| TYPE | Bone Levers |
| Aesculap | FK185 |
| KLS Martin | 23-541-08 |
| AUTHORS | HOHMANN |
| Medicon | 77.11.10 |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|-----------------------------|
| UNIT | pcs. |
| TYPE | Bone Levers |
| Aesculap | FK185 |
| KLS Martin | 23-541-08 |
| AUTHORS | HOHMANN |

Medicon

77.11.10