

HOHMANN BONE LEVER

43mm; narrow tip

24.0cm (9.6")


SKU: 190-572-43

Categories: [Open Surgery](#), [Orthopaedics](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bone Levers
Aesculap	FK188
KLS Martin	23-543-43
AUTHORS	HOHMANN



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bone Levers
Aesculap	FK188
KLS Martin	23-543-43
AUTHORS	HOHMANN