

## HOHMANN BONE LEVER

24mm; wide, round end; 1 whole handle

27.0cm (10.8")


**SKU:** 190-573-24

**Categories:** [Open Surgery](#), [Orthopaedics](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bone Levers</a>
<b>Aesculap</b>	FK176R
<b>KLS Martin</b>	23-545-22
<b>AUTHORS</b>	<a href="#">HOHMANN</a>
<b>Medicon</b>	77.10.95



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bone Levers</a>
<b>Aesculap</b>	FK176R
<b>KLS Martin</b>	23-545-22
<b>AUTHORS</b>	<a href="#">HOHMANN</a>

**Medicon**

77.10.95