

HOHMANN BONE LEVER

24mm; wide, round end; 1 whole handle

27.0cm (10.8")


SKU: 190-573-24

Categories: [Open Surgery](#), [Orthopaedics](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bone Levers
Aesculap	FK176R
KLS Martin	23-545-22
AUTHORS	HOHMANN
Medicon	77.10.95



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bone Levers
Aesculap	FK176R
KLS Martin	23-545-22
AUTHORS	HOHMANN

Medicon

77.10.95