

HOHMANN BONE LEVER

35mm; wide tip

29.0cm (11.6")


SKU: 190-573-35

Categories: [Open Surgery](#), [Orthopaedics](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bone Levers
Aesculap	FK173
KLS Martin	23-545-35
AUTHORS	HOHMANN
Medicon	77.11.34



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bone Levers
Aesculap	FK173
KLS Martin	23-545-35
AUTHORS	HOHMANN

Medicon

77.11.34