

HOHMANN BONE LEVER

43mm; narrow tip; 1 whole handle

24.0cm (9.6")


SKU: 190-573-43

Categories: [Open Surgery](#), [Orthopaedics](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bone Levers
Aesculap	FK174
KLS Martin	23-541-43
AUTHORS	HOHMANN
Medicon	77.11.42



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bone Levers
Aesculap	FK174
KLS Martin	23-541-43
AUTHORS	HOHMANN

Medicon

77.11.42