

HOLLOW OBTURATOR FOR NEPHROSCOPY SHEATH

23 Fr.; 630-423-12; 15°; FL 220mm

drilling Ø 1.4mm

fits Compact-Nephroscopes 601-135-45 + 601-135-90 Colour Code: white

SKU: 630-423-17

Categories: [Endoscopy](#), [Urology](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
RANGE [Nephroscopy](#)
TYPE [Obtulators](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
RANGE [Nephroscopy](#)
TYPE [Obtulators](#)