

## HOLLOW OBTURATOR FOR NEPHROSCOPY SHEATH

23 Fr.; 630-423-12; 15°; FL 220mm

drilling Ø 1.4mm

fits Compact-Nephroscopes 601-135-45 + 601-135-90 Colour Code: white

**SKU:** 630-423-17

**Categories:** [Endoscopy](#), [Urology](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)  
**RANGE** [Nephroscopy](#)  
**TYPE** [Obturers](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)  
**RANGE** [Nephroscopy](#)  
**TYPE** [Obturers](#)