

HOOK ELECTRODE; 9 FR.

retrograd angled

for Ø 1.9mm endoscope; 0°



SKU: 630-512-70

Categories: [Endoscopy](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Pediatric Resectoscopy
TYPE	Electrodes

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Pediatric Resectoscopy
TYPE	Electrodes