

HTT FIBER LIGHT GUIDE CABLE

for Xenon light source; w/o connectors

Ø 2.5 x 3000 mm Length

tube colour yellow

HTT
HIGH TEMPERATURE TECHNOLOGY LINE



SKU: 540-812-30YW

Categories: [Arthroscopy](#), [Endoscopy](#), [Equipment](#), [GUBBINI/Gyn](#), [Gynecology](#), [Laparoscopy](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Light Guide Cables](#)

TYPE [HTT \(High Temperature Technology\)](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Light Guide Cables](#)

TYPE [HTT \(High Temperature Technology\)](#)