


HTT FIBER LIGHT GUIDE CABLE

for Xenon light source; w/o connectors

Ø 3.5 x 1800 mm Length

Ø 3.5 x 1800 mm Length, Tube colour green

HTT
HIGH TEMPERATURE TECHNOLOGY LINE



SKU: 540-813-18GN

Categories: [Arthroscopy](#), [Endoscopy](#), [Equipment](#), [GUBBINI/Gyn](#), [Gynecology](#), [Laparoscopy](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Light Guide Cables](#)

TYPE [HTT \(High Temperature Technology\)](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Light Guide Cables](#)

TYPE [HTT \(High Temperature Technology\)](#)