


## HTT FIBER LIGHT GUIDE CABLE

for Xenon light source; w/o connectors

Ø 3.5 x 1800 mm Length

Tube colour lilac

**HTT**  
HIGH TEMPERATURE TECHNOLOGY LINE



**SKU:** 540-813-18LL

**Categories:** [Arthroscopy](#), [Endoscopy](#), [Equipment](#), [GUBBINI/Gyn](#), [Gynecology](#), [Laparoscopy](#), [Urology](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Light Guide Cables</a>
<b>TYPE</b>	<a href="#">HTT (High Temperature Technology)</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Light Guide Cables</a>
<b>TYPE</b>	<a href="#">HTT (High Temperature Technology)</a>