

## HTT FIBER LIGHT GUIDE CABLE

for Xenon light source; w/o connectors

Ø 3.5 x 2300 mm Length

tube colour yellow

**HTT**  
HIGH TEMPERATURE TECHNOLOGY LINE



**SKU:** 540-813-23YW

**Categories:** [Arthroscopy](#), [Endoscopy](#), [Equipment](#), [GUBBINI/Gyn](#), [Gynecology](#), [Laparoscopy](#), [Urology](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Light Guide Cables](#)

**TYPE** [HTT \(High Temperature Technology\)](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Light Guide Cables](#)

**TYPE** [HTT \(High Temperature Technology\)](#)