


## HTT FIBER LIGHT GUIDE CABLE

for Xenon light source; w/o connectors

Ø 4.9 x 5000 mm Length

tube colour lilac

**HTT**  
HIGH TEMPERATURE TECHNOLOGY LINE



**SKU:** 540-815-50LL

**Categories:** [Arthroscopy](#), [Endoscopy](#), [Equipment](#), [GUBBINI/Gyn](#), [Gynecology](#), [Laparoscopy](#), [Urology](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Light Guide Cables</a>
<b>TYPE</b>	<a href="#">HTT (High Temperature Technology)</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Light Guide Cables</a>
<b>TYPE</b>	<a href="#">HTT (High Temperature Technology)</a>