

HTT FIBER LIGHT GUIDE CABLE

for Xenon light source; w/o connectors

Ø 2.5 x 2300 mm Length

tube colour lilac

HTT
HIGH TEMPERATURE TECHNOLOGY LINE



SKU: 540-812-23LL

Categories: [Arthroscopy](#), [Endoscopy](#), [Equipment](#), [GUBBINI/Gyn](#), [Gynecology](#), [Laparoscopy](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Light Guide Cables
TYPE	HTT (High Temperature Technology)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Light Guide Cables
TYPE	HTT (High Temperature Technology)