


HTT FIBER Y LIGHT GUIDE CABLE

for Xenon light source; w/o connectors

Y-Typ 1800mm; 3.5/4.8; 2x Proj / 1x Endo

tube colour lilac

HTT
HIGH TEMPERATURE TECHNOLOGY LINE



SKU: 540-817-18LL

Categories: [Arthroscopy](#), [Endoscopy](#), [Equipment](#), [GUBBINI/Gyn](#), [Gynecology](#), [Laparoscopy](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|--------------|---|
| UNIT | pcs. |
| RANGE | Light Guide Cables |
| TYPE | HTT (High Temperature Technology) |

ZUSÄTZLICHE INFORMATIONEN

| | |
|--------------|---|
| UNIT | pcs. |
| RANGE | Light Guide Cables |
| TYPE | HTT (High Temperature Technology) |