

## HTT FIBER Y LIGHT GUIDE CABLE

for Xenon light source; w/o connectors

Y-Typ 2300mm; 3.5/4.8; 2x Proj / 1x Endo

tube colour grey

HTT

HIGH TEMPERATURE TECHNOLOGY LINE



**SKU:** 540-817-23

**Categories:** [Arthroscopy](#), [Endoscopy](#), [Equipment](#), [GUBBINI/Gyn](#), [Gynecology](#), [Laparoscopy](#), [Urology](#)

### ZUSÄTZLICHE INFORMATIONEN

|              |                                                   |
|--------------|---------------------------------------------------|
| <b>UNIT</b>  | <a href="#">pcs.</a>                              |
| <b>RANGE</b> | <a href="#">Light Guide Cables</a>                |
| <b>TYPE</b>  | <a href="#">HTT (High Temperature Technology)</a> |

### ZUSÄTZLICHE INFORMATIONEN

|              |                                                   |
|--------------|---------------------------------------------------|
| <b>UNIT</b>  | <a href="#">pcs.</a>                              |
| <b>RANGE</b> | <a href="#">Light Guide Cables</a>                |
| <b>TYPE</b>  | <a href="#">HTT (High Temperature Technology)</a> |