

## HYSTEROSCOPE HD+

Ø 4.0mm; 70°; 300 mm

autoclavable


**SKU:** 601-114-70HD

**Categories:** [Endoscopy](#), [Equipment](#), [Gynecology](#)



### ZUSÄTZLICHE INFORMATIONEN

|              |                            |
|--------------|----------------------------|
| <b>UNIT</b>  | <a href="#">pcs.</a>       |
| <b>RANGE</b> | <a href="#">Telescopes</a> |
| <b>SCOPE</b> | <a href="#">Gynecology</a> |



### ZUSÄTZLICHE INFORMATIONEN

|              |                            |
|--------------|----------------------------|
| <b>UNIT</b>  | <a href="#">pcs.</a>       |
| <b>RANGE</b> | <a href="#">Telescopes</a> |
| <b>SCOPE</b> | <a href="#">Gynecology</a> |