

## HYSTEROSCOPE HD+2

Ø 4.0mm; 12°; 300 mm

autoclavable

**SKU:** 601-114-12HD2

**Categories:** [Endoscopy](#), [Equipment](#), [Gynecology](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Telescopes</a>
<b>SCOPE</b>	<a href="#">Gynecology</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Telescopes</a>
<b>SCOPE</b>	<a href="#">Gynecology</a>