

## HYSTEROSCOPE HD+2

Ø 4.0mm; 70°; 300 mm

autoclavable

**SKU:** 601-114-70HD2

**Categories:** [Endoscopy](#), [Equipment](#), [Gynecology](#)




### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Telescopes](#)

**SCOPE** [Gynecology](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Telescopes](#)

**SCOPE** [Gynecology](#)