

HYSTEROSCOPY DIAGNOSTIC SHEATH

0°; 1 turnable stopcock

Ø 2.6 mm; 275 mm

SKU: 640-010-27

Categories: [Endoscopy](#), [Gynecology](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Hysteroscopy
TYPE	Diagnostic Sheaths

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Hysteroscopy
TYPE	Diagnostic Sheaths