

HYSTEROSCOPY DIAGNOSTIC SHEATH

0°; 1 turnable stopcock

Ø 4.0 mm; 270 mm; white



SKU: 640-010-30

Categories: [Endoscopy](#), [Gynecology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Hysteroscopy
TYPE	Diagnostic Sheaths

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Hysteroscopy
TYPE	Diagnostic Sheaths