

HYSTEROSCOPY DIAGNOSTIC SHEATH

30°; 1 turnable stopcock

Ø 4.0 mm; 270 mm; white

for telescope 2,9mm/30°/303mm



SKU: 640-010-33

Categories: [Endoscopy](#), [Gynecology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Hysteroscopy
TYPE	Diagnostic Sheaths

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Hysteroscopy
TYPE	Diagnostic Sheaths