

## HYSTEROSCOPY DIAGNOSTIC SHEATH

0°; 1 turnable stopcock

Ø 5.0 mm; 270 mm; blue



**SKU:** 640-010-40

**Categories:** [Endoscopy](#), [Gynecology](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Hysteroscopy</a>
<b>TYPE</b>	<a href="#">Diagnostic Sheaths</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Hysteroscopy</a>
<b>TYPE</b>	<a href="#">Diagnostic Sheaths</a>