

## HYSTEROSCOPY DIAGNOSTIC SHEATH

0°; 1 turnable stopcock

Ø 5.0 mm; 270 mm; blue

**SKU:** 640-010-40

**Categories:** [Endoscopy](#), [Gynecology](#)



### ZUSÄTZLICHE INFORMATIONEN

|              |                                    |
|--------------|------------------------------------|
| <b>UNIT</b>  | <a href="#">pcs.</a>               |
| <b>RANGE</b> | <a href="#">Hysteroscopy</a>       |
| <b>TYPE</b>  | <a href="#">Diagnostic Sheaths</a> |

### ZUSÄTZLICHE INFORMATIONEN

|              |                                    |
|--------------|------------------------------------|
| <b>UNIT</b>  | <a href="#">pcs.</a>               |
| <b>RANGE</b> | <a href="#">Hysteroscopy</a>       |
| <b>TYPE</b>  | <a href="#">Diagnostic Sheaths</a> |