

HYSTEROSCOPY DIAGNOSTIC SHEATH

30°; 1 turnable stopcock

Ø 5.0 mm; 270 mm; blue



SKU: 640-010-43

Categories: [Endoscopy](#), [Gynecology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Hysteroscopy
TYPE	Diagnostic Sheaths

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Hysteroscopy
TYPE	Diagnostic Sheaths