

HYSTEROSCOPY DIAGNOSTIC SHEATH

30°; 1 turnable stopcock

Ø 5.0 mm; 270 mm; blue



SKU: 640-010-43

Categories: [Endoscopy](#), [Gynecology](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|--------------|------------------------------------|
| UNIT | pcs. |
| RANGE | Hysteroscopy |
| TYPE | Diagnostic Sheaths |

ZUSÄTZLICHE INFORMATIONEN

| | |
|--------------|------------------------------------|
| UNIT | pcs. |
| RANGE | Hysteroscopy |
| TYPE | Diagnostic Sheaths |