

HYSTEROSCOPY DIAGNOSTIC SHEATH

0°; 1 turnable stopcock

Ø 2.9 mm; 228 mm



SKU: 640-010-20

Categories: [Endoscopy](#), [Gynecology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Hysteroscopy
TYPE	Diagnostic Sheaths

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Hysteroscopy
TYPE	Diagnostic Sheaths