

## INNER TUBE ONLY FOR BIPOLAR

coagulation (ex.628-201-53)

forceps 628-201-53 and 628-205-53; 340mm

**SKU:** 628-204-53

**Categories:** [Endoscopy](#), [Laparoscopy](#)



## ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">HF-Instruments, bipolar</a>
<b>TYPE</b>	<a href="#">Sheaths</a>

## ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">HF-Instruments, bipolar</a>
<b>TYPE</b>	<a href="#">Sheaths</a>