

JANET WOUND AND BLADDER SYRINGE

with two exchangeable tips; 50 ml

SKU: 087-070-05

Categories: [Open Surgery](#), [Puncture / Special Needles](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Syringes
Aesculap	SH405
KLS Martin	18-930-05
AUTHORS	JANET
Medicon	09.78.81

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Syringes
Aesculap	SH405
KLS Martin	18-930-05
AUTHORS	JANET
Medicon	09.78.81