

## JANET WOUND AND BLADDER SYRINGE

with two exchangeable tips; 50 ml

**SKU:** 087-070-05

**Categories:** [Open Surgery](#), [Puncture / Special Needles](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Syringes</a>
<b>Aesculap</b>	SH405
<b>KLS Martin</b>	18-930-05
<b>AUTHORS</b>	<a href="#">JANET</a>
<b>Medicon</b>	09.78.81

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Syringes</a>
<b>Aesculap</b>	SH405
<b>KLS Martin</b>	18-930-05
<b>AUTHORS</b>	<a href="#">JANET</a>
<b>Medicon</b>	09.78.81