

JANSEN EAR FORCEPS

bayonet shaped; serrated

20.0cm (8"); width of jaw 2.6 +/-0.2



SKU: 020-760-20

Categories: [Forceps](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Ear and Nose Dressing Forceps
KLS Martin	36-395-20
AUTHORS	JANSEN
SPECIFICATIONS	bajonett , serrated



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Ear and Nose Dressing Forceps
KLS Martin	36-395-20
AUTHORS	JANSEN
SPECIFICATIONS	bajonett , serrated