

## JANSEN EAR FORCEPS

bayonet shaped; serrated

20.0cm (8"); width of jaw 2.6 +/-0.2



**SKU:** 020-760-20

**Categories:** [Forceps](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Ear and Nose Dressing Forceps</a>
<b>KLS Martin</b>	36-395-20
<b>AUTHORS</b>	<a href="#">JANSEN</a>
<b>SPECIFICATIONS</b>	<a href="#">bajonett</a> , <a href="#">serrated</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Ear and Nose Dressing Forceps</a>
<b>KLS Martin</b>	36-395-20
<b>AUTHORS</b>	<a href="#">JANSEN</a>
<b>SPECIFICATIONS</b>	<a href="#">bajonett</a> , <a href="#">serrated</a>