

JANSEN EAR FORCEPS

bayonet shaped; serrated

16.0cm (6.4"); width of jaw 2.0 +/-0.2



SKU: 020-760-16

Categories: [Forceps](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Ear and Nose Dressing Forceps
Aesculap	BD880
KLS Martin	36-395-16
AUTHORS	JANSEN
Medicon	07.30.16
SPECIFICATIONS	bajonett , serrated



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Ear and Nose Dressing Forceps
Aesculap	BD880
KLS Martin	36-395-16

AUTHORS

[JANSEN](#)

Medicon

07.30.16

SPECIFICATIONS

[bajonett](#), [serrated](#)